

Safeguarding the water quality

The simple and quick solution: complete filter

Pre-filter (optional)
Complete filter with 1,5 µm retention rate



Main filter
Complete filter with 0,2 µm sterile filtration membrane



Endotoxin retention
Complete filter 0,2 µm with endotoxin removal




Secure cleaning and disinfection device




The alternative for complex systems

Stainless steel housing with particle filter 1,5 µm



Step 2 Variant A

Stainless steel housing with sterile filter 0,2 µm



Step 2 Variant B

Stainless steel housing with sterile filter 0,2 µm with endotoxin removal

